

# Strategies for the Elimination of Mycoplasma Contamination from Laboratory Incubators

## Introduction:

Mycoplasma contamination poses a significant challenge in cell culture laboratories, compromising research integrity and leading to inaccurate experimental results. Estimates for the incidence of Mycoplasma contamination vary from 10% to 36% of cell lines used in laboratory procedures. Laboratory incubators, being a critical environment for cell growth, are prone to Mycoplasma contamination. This paper explores effective strategies to eliminate Mycoplasma contamination from laboratory incubators, ensuring reliable and consistent cell culture experiments.

# I. Cleaning and Decontamination Procedures:

To eliminate *Mycoplasma* contamination, thorough cleaning and decontamination protocols are essential. The following steps can be employed:

- 1. Regular Cleaning: Implement a stringent cleaning regimen for laboratory incubators, involving the removal of all cell culture vessels, shelves, and trays. Clean surfaces using appropriate disinfectants, such as 70% ethanol or other proven disinfectants, to effectively eliminate mycoplasma.
- 2. Disinfection: After cleaning, disinfect the incubator using a sterilizing agent, such as nebulized hydrogen peroxide (nHP), Paracetic acid (PA) or ultraviolet (UV) light. These methods are highly effective in eradicating *Mycoplasma* and other contaminants.
- 3. Periodic Maintenance: Establish a maintenance schedule for incubators to ensure proper functioning and prevent *Mycoplasma* contamination. This includes regular replacement of filters, inspection of seals, and calibration of temperature and humidity controls.

#### II. Isolation and Quarantine:

Preventing *Mycoplasma* contamination requires isolating potentially contaminated cell lines and practicing strict quarantine measures:

1. Isolation: Identify and segregate *Mycoplasma* -contaminated cell lines from healthy ones. Store contaminated cultures separately, preferably in sealed containers, to prevent cross-contamination.



2. Quarantine: Newly obtained or suspected cell lines should be quarantined for *Mycoplasma* testing before integration into the laboratory workflow. Maintain a dedicated space for quarantine and perform regular testing to ensure early detection and prompt elimination of contamination sources.

# III. Regular Testing:

Regular and rigorous testing protocols are indispensable in identifying mycoplasma contamination and taking appropriate measures:

- 1. Testing Methods: Employ validated *Mycoplasma* detection methods, such as polymerase chain reaction (PCR) or DNA staining assays. These tests are highly sensitive and specific, allowing for accurate detection of *Mycoplasma* contamination.
- 2. Frequency: Establish a periodic testing schedule to ensure the timely identification of *Mycoplasma* contamination. Regular testing, at least once a month, is recommended, especially for high-throughput laboratories with extensive cell culture activities.
- 3. External Testing: Periodically utilize external services or contract laboratories for mycoplasma testing to obtain an independent and unbiased assessment of incubator and cell culture contamination.

## Conclusion:

Mycoplasma contamination represents a persistent challenge in laboratory incubators and can compromise the validity of cell culture experiments. By implementing robust cleaning and decontamination procedures, ensuring isolation and quarantine practices, and implementing regular Mycoplasma testing, laboratories can significantly reduce the risk of contamination. Maintaining a vigilant approach towards Mycoplasma elimination is crucial for generating reliable and reproducible experimental results, thus promoting the advancement of scientific research.

### References

Drexler H, Uphof CC (2002). In Vitro Cell. Biol. - Animal 38:79-85

Olarerin-George AO, Hogenesch JB (2015). Nucleic Acids Research 43(5): 2535-2542

Chat-GPT 4, OpenAl