



Order Number

75400-714

Capacity

0.52 cu.ft. (14.8 liters)

Interior Dimensions (W x D x H)

12.5" x 6.1" x 6.0"
318 x 154 x 152 mm

Exterior Dimensions (W x D x H)

30.4" x 24.2" x 12.2"
772 x 615 x 310 mm

Dry Weight*

80.8 lb (36.7 kg)

Net Weight**

125.6 lb (57.0 kg)

Shipping Weight

167.0 lb (76.0 kg)

Interior Construction

Aluminum

Exterior Construction

300 Stainless Steel/CR

Wall Clearance, Each Side

4" (102 mm)

Wall Clearance, Rear

4" (102 mm)

Tank Size

14 liters

Tank Fill Capacity

10 liters

Electrical

Electrical Supply /Hz, AC, 1Ø

50/60 Hz,
100-240 AC,
550 Watts

Plug Supplied

NEMA5-15P, EU1-16P

Temperature Specifications

Temperature Range

2°C to 70°C

Temperature Uniformity

±0.5°C at 4°C
±0.5°C at 37°C
±1.3°C at 70°C

Temperature Stability - At 37°C

0.1

*Bath empty with no beads loaded.

**Bath filled with beads.

We reserve the right to alter specifications at any time

