



co-culture duet system

technical specifications sheet



general

compatible products	stratus and other standard plate readers that can accept a 15.0mm tall microplate (~17.5mm with a lid)
capacity	up to 18 duets
duet	2 reservoirs, 1 membrane
dimensions	frame: 27.6 x 85.05 x 15.0mm (LWH) reservoir, internal: 8 x 9 x 12mm (LWH) duet: 24 x 15 x 14mm (LWH)
weight	frame: 89g duets: 3g each
reservoir volume	800-1000 μ L

membrane

material	polycarbonate
treatment	PVP (hydrophilic)
pore size	0.2 μ m
surface area	50.3mm ²
diameter	8mm
membrane thickness	25 μ m

All cerillo products are for research use only and are not intended for use in diagnostic procedures.

To learn more about Cerillo and our products, visit us at www.cerillo.bio or by email at info@cerillo.bio