



canopy

technical specifications sheet

The canopy is RoHS, REACH and FCC compliant, CE certified, made in the USA, and has a two-year warranty.



general

compatible products	stratus microplate reader
dimensions	113.2 X 66.1 X 47.8mm (LWH)
weight	0.17 kg
power supply	power = 3A, voltage = 5.1V
input	USB-C USB-C 5.1 V DC, 3A

communication

wired communication	USB 2.0
wireless communication range	15 m

cerillo labrador premium software

operating system	Windows: 10, 11 MacOS: 12 (Monterey)
data export format	CSV

All cerillo products are for research use only and are not intended for use in diagnostic procedures.

To learn more about cerillo and our products , visit us at www.cerillo.bio or email us at info@cerillo.bio