



**KENTRON**  
MICROBIOLOGY

*"The A-Box, a square anaerobic jar, that makes your anaerobic cultures much simpler."*



*Looking for a  
square anaerobic jar?  
Try the A-Box*



**KENTRON**  
MICROBIOLOGY

*"Make your lab-life easier  
with the A-box."*

## Looking for a square anaerobic jar?

If you are looking for a square anaerobic jar that fits micro-titer plates for your microbiology cultures under anaerobic conditions, then look no further.

We proudly present to you the A-Box. This unique product is designed to use with microtiter plates (deep wells) for your cultures under anaerobic conditions, large testtubes or petri-dishes. The A-box can be used with multiple anaerobic systems, Anaerobic Chambers and GasPak methods.

Finally the moment we have all been waiting for, it's here. The A-Box, a square anaerobic jar, that makes your anaerobic cultures much simpler.



## Technical specifications:

External dimensions: 18.0 x 18.4 x 18.6 cm

Internal dimensions: 14.0 x 17.0 x 13.1 cm

Housing: PMMA Lid: stainless steel, lid contains 2 clips for catalysts. The A-box will be supplied with 1 rack for microtiterplates (test tube racks optional).



**Kentron Microbiology B.V.**

Logistiekweg 20-27

7007 CJ Doetinchem The Netherlands

Tel: +31-314-380527

Fax: +31-314-382899

Email: [office@kentron-microbiology.nl](mailto:office@kentron-microbiology.nl)

Website: [www.kentron-microbiology.nl](http://www.kentron-microbiology.nl)